


Chiba

Chiba tip, disposable needle
Aiguille à usage unique avec pointe Chiba
Aguja desechable con punta Chiba



- The Chiba model disposable needle is designed for cytological aspiration, amniocentesis, and local anaesthesia. It features a cannula with centimetre markings, an echogenic treatment ensuring precise positioning under ultrasound guidance, and adheres to an internationally recognized colour code for diameter identification.
- L'aiguille modèle Chiba à usage unique est conçue pour l'aspiration cytologique, l'amniocentèse et l'anesthésie locale. La canule est dotée d'un marquage centimétrique, d'un traitement échogène permettant un positionnement précis sous guidage échographique, ainsi que d'un code couleur internationalement reconnu pour l'identification du diamètre.
- La aguja desechable modelo Chiba está diseñada para aspiración citológica, amniocentesis y anestesia local. La cánula tiene un marcado centimétrico, un tratamiento ecogénico que permite un correcto posicionamiento bajo guía ecográfica y un código de colores reconocido internacionalmente para la identificación del diámetro.

Available sizes
Tailles disponibles
Medidas disponibles

	60 mm	80 mm	100 mm	120 mm	160 mm	180 mm	200 mm	250 mm	300 mm
18G	CH1806	CH1808	CH1810	CH1812	CH1816	CH1818	CH1820	CH1825	CH1830
20G	CH2006	CH2008	CH2010	CH2012	CH2016	CH2018	CH2020	CH2025	CH2030
21G	CH2106	CH2108	CH2110	CH2112	CH2116	CH2118	CH2120	CH2125	CH2130
22G	CH2206	CH2208	CH2210	CH2212	CH2216	CH2218	CH2220	CH2225	CH2230

Available in boxes of 20 pcs box (supplied in single packaging) - EtO sterilization - CE 0051 marked
 Disponible en boîtes de 20 pièces (emballé sous sachet individuel) - stérilisation EtO - marqué CE 0051
 Disponible en cajas de 20 piezas (se suministra en un solo paquete) - esterilización por EtO - marcado CE 0051